

HANDLING TUBE

Implementing animal welfare while handling mice with a tube.

“Research has shown that picking up laboratory mice by the tail induces aversion and high anxiety in these animals, and generally should be avoided. Instead, wherever possible, mice should be picked up by non-aversive methods such as tunnel handling and cupping, which promote positive responses to human contact “.

Source : <https://www.nc3rs.org.uk/mouse-handling-tutorial>



Reduce anxiety in mice

Mice show more reliable behavioural and physiological responses

Mice only need brief experience of tunnel handling to habituate

Tunnel handling becomes as easy as tail handling within short time

SERLAB UK Ltd

Lincoln House - ccfgb

300 High Holborn

London WC1V 7JH

Phone : +44 (0) 207 0923436

Email : sales@serlab.co.uk

www.serlab.co.uk

DESCRIPTION

Material : clear high quality polycarbonate

Dimensions : L 130mm - diameter 52mm - thickness 1.5mm

Autoclavable - cage washer safe

Carton box of 54 units