

LAB ASPEN G7

DESCRIPTION

Consistent cross cut chips obtained through the mechanical chopping of non chemically treated Aspen logs. Sanitized treated at high temperature through a drum dryer.



PHYSICAL PROPERTIES

Origin : France
Raw material : non chemically treated aspen logs
pH : 5,0 to 5,6
Density : 0.175 to 0.190
Moisture content : 13 % maximum
Particule size : 80% min. between 2 and 4 mm
Absorption : 260 %

SERLAB UK Ltd

Lincoln House - ccfgb

300 High Holborn

London WC1V 7JH

Phone : +44 (0) 207 0923436

Email : sales@serlab.co.uk

www.serlab.co.uk

PACKAGING

10 kg autoclavable polypropylene sack
Every bag labelled with a batch number
39 sacks per shrink wrapped pallet

GLP Certificate
Irradiation : optional