


SIZZLE-dri

DESCRIPTION

Nesting material made of paper strands.
Available in brown or white.



PHYSICAL PROPERTIES

Material : Kraft Paper.
Colour : Brown or white.
Density: 65 kg/m³
Moisture content : 8 à 12 %

SERLAB UK Ltd

Lincoln House - ccfgb

300 High Holborn

London WC1V 7JH

Phone : +44 (0) 207 0923436

Email : sales@serlab.co.uk

www.serlab.co.uk

PACKAGING

Autoclavable 5 kg bag in a cardboard box.
Pallet size : 24 boxes - 80 x 120.

Contaminant screened.
Irradiation : optional